

Ophthalmic unit



■ Features

- Fashion dedign
- Table with rotate and motorized up/down function




Table motorized up/down



Table rotate



Rotatable tray table for install the Auto Ref-keratometer



Motorized up/down chair and phoropter arm



Specification

Specification		
Table	Dimensions	1040(W)*1030(D)*1830(H)mm
	Table size	1000(W)*480(D)mm
	Lift range for phoropter arm	420mm
	Power supply	AC 220V/55Hz,110V/60Hz,120W
	Net weight	150.0Kg
Chair	Lift range for chair	480~630mm
	Net weight	45.0Kg